ASATEX®

Art. **PL-PG**

Sizes **9 - 10**

€ **0598**

PPE - CAT. III - High risks

EN ISO 21420:2020

EN 388:2016+A1:2018

EN ISO 374-1:2016+A1:2018 / Type A

EN ISO 374-4:2019

EN ISO 374-5:2016

Dexterity: 5

EN 388:2016
+A1:2018

4 1 2 1 X

EN ISO 374-1:2016

TYPE A

J K L M P S T

EN ISO 374-5:2016

VIRUS

MM.YYYY

(DE) Herstellerinformation

(EN) Manufacturer's Information

(ES) Información del fabricante

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Manufacturer's Information (EN)

according to Regulation (EU) 2016/425, Annex II, Section 1.4. (published in the Official Journal of the European Union). Please read carefully before using! You are required to enclose this information leaflet when passing on the personal protective equipment (PPE), or to present it personally to the recipient. You may therefore reproduce this leaflet at your own discretion.

Declaration of Conformity: These gloves are classified as personal protective equipment (PPE). The CE mark confirms that the product satisfies the applicable requirements of Regulation (EU) 2016/425. For the full Declaration of Conformity, please visit: www.asatex.eu/konf

A. Explanation and numbers of the standards whose requirements the gloves satisfy: Standards retrieved from the Official Journal of the European Union. Available from Beuth Verlag GmbH, 10787 Berlin. www.beuth.de.

EN ISO 21420:2020 - General requirements and test methods for gloves. The verification of dexterity, performed according to the EN ISO 21420:2020 see page 1.

EN 388:2016+A1:2018 - Protective gloves against mechanical risks

EN 388:2016 +A1:2018	Test criteria	Rating
	A = Abrasion resistance	0 - 4
	B = Cut resistance (Coupe test)	0 - 5
	C = Tear resistance	0 - 4
	D = Puncture resistance	0 - 4
	E = Cut resistance (TDM) according to EN ISO 13997:1999	A - F
	F = Impact protection test according to EN 13594:2015	P

The higher the test number, the better the test performance. X means 'not tested'; P means 'passed'. The result relates to palm of the hand.

EN ISO 374-1:2016+A1:2018 - Protective gloves against micro-organisms and dangerous chemicals.

EN ISO 374-1:2016+A1:2018 / Type A - Part 1: Terminology and performance requirements for chemical risks.

EN ISO 374-1:2016 TYPE A	EN ISO 374-1:2016 TYPE B	EN ISO 374-1:2016 TYPE C
Penetration resistance (EN ISO 374-2) Breakthrough time 30 min for minimum 3 products of the list (EN 16523-1).	Penetration resistance (EN ISO 374-2) Breakthrough time 30 min for minimum 3 products of the list (EN 16523-1).	Penetration resistance (EN ISO 374-2) Breakthrough time 10 min for minimum 1 product of the list (EN 16523-1).

Level	1	2	3	4	5	6
Breakthrough time (min.)	>10	>30	>60	>120	>240	>480

The permeation resistance has been assessed under laboratory conditions and relates only to tested specimen.

EN ISO 374-2:2019 Part 2: Determination of resistance to penetration. Penetration test.
EN ISO 374-4:2019 Part 4: Determination of resistance to degradation by chemicals. Degradation test.
Degradation results indicate the changes in puncture resistance of the gloves after exposure to the challenged Chemical.

EN ISO 374-5:2016 Protective gloves against chemicals and dangerous micro-organisms.
Part 5: Terminology and performance requirements for risks by micro-organisms.
The penetration resistance has been assessed under laboratory condition and relates only to the tested specimen.

EN ISO 374-5:2016	EN ISO 374-5:2016
 VIRUS	Gloves who offer protection against bacteria, fungi and micro-organisms.
	Gloves who offer protection against bacteria and fungi.

PERMEATION DATA			
Code letter	Test chemical	CAS-RN	Class
A	Methanol	67-56-1	Primary alcohol
B	Acetone	67-64-1	Ketone
C	Acetonitrile	75-05-8	Nitrile
D	Dichloromethane	75-09-2	Chlorinated hydrocarbon
E	Carbon sulphide	75-15-0	Sulphur-containing organic compound
F	Toluene	108-88-3	Aromatic hydrocarbon
G	Diethylamine	109-89-7	Amine
H	Tetrahydrofuran	109-99-9	Heterocyclic and ether compounds
I	Ethyl acetate	141-78-6	Ester
J	n-heptane	142-82-5	Aliphatic hydrocarbons
K	Sodium hydroxide 40%	1310-73-2	Inorganic alkali
L	Sulphuric acid 96%	7664-93-9	Inorganic acid, oxidizing
M	Nitric acid 65%	7697-37-2	Inorganic acid, oxidizing
N	Acetic acid 99%	64-19-7	Organic acid
O	Ammonia water 25%	1336-21-6	Organic alkali
P	Hydrogen peroxide 30%	7722-84-1	Peroxide
S	Hydrofluoric acid 40%	7664-39-3	Inorganic acid
T	Formaldehyde 37%	50-00-0	Aldehyde

WARNINGS:

- This information does not provide any details on the actual duration of protection at the workplace; it also does not distinguish between blends and pure chemicals.
- Resistance to chemicals was assessed using samples taken only from the palm and tested under laboratory conditions (apart from the glove measures 400 mm, in which case the cu is also tested); the stated resistance refers only to the tested chemicals. Resistance may differ if the chemical is present in a blend.
- Users are recommended to check whether the glove is suitable for its intended application, as the conditions at the workplace may differ from those during type testing, depending on the temperature, abrasion and degradation.
- Protective gloves that have already been used may provide less resistance to dangerous chemicals due to changes in their physical properties. The actual service life may be reduced significantly due to degradation, movement, stringing, abrasion and suchlike, caused by contact with chemicals. Degradation may be the most significant factor in regard to aggressive chemicals; this must be duly considered in the selection of protective gloves against chemicals.
- The gloves must always be checked for imperfections before use. Glove can be cleaned by wiping with damp cloth.

B. Purpose, applications and risk assessment:

These gloves are intended for protection against medium risks. These products satisfy the requirements of the quoted standards. Please note that the actual conditions of use cannot be simulated and that the decision on the products suitability for its intended purpose therefore lies exclusively with the user. The manufacturer is not responsible for improper use. Hence, an assessment of the residual risk should be performed before use in order to determine whether this product is suitable for its intended purpose. Please guide by the pictograms and the indicated performance levels. Due to its resistance to oils and greases, PVC is ideal for use in refineries, when working with fuels or when handling oils and greases.

C. Cleaning, care and disinfecting - Not advised for cleaning

Precautionary measures during use

- Do not use these gloves to protect against serrated edges or blades.
- If gloves must be used in a hot environment, make certain that they satisfy the requirements of EN 407 and that they were tested as specified therein.
- The gloves should not be used in the vicinity of moving machine parts due to the risk of entrapment.
- Check the gloves carefully before use to make certain there are no defects or imperfections.
- If gloves fulfil the requirements for puncture resistance in accordance with EN 388:2016, this does not automatically mean that they also offer protection against puncture by sharp objects such as hypodermic needles. The puncture protection does not offer protection against hypodermic needles.
- Don't use gloves on the inside as they may lead to skin irritation due to the dirty surface. Used gloves may have a lower resistance against dangerous chemicals. Don't reuse contaminated gloves.

D. Storage and aging - Keep in a cool, dry place; do not expose to direct sunlight; keep away from any ignition sources; store in the original packaging if possible. The mechanical properties of the gloves will not change for a period of up to three years from the manufacturing date, provided they are stored as recommended. The manufacturing date (month/year) is stated on the glove.

E. Disposal - Discard gloves according to current regulations.
The manufacturer was examined under the system of ensuring EU Quality of production by means of monitoring EU Regulation 2016/425 Module D SGS Fimko Oy P.O. Box 30 (Särkinientie 3) 00211 Helsinki Finnland - Notified Body number - 0598

For additional technical information, please visit: www.asatex.eu

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